

CCDC65 antibody



Catalog Number: 108984

Product name

CCDC65 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CCDC65 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 57 kDa

Preparation

This antibody was obtained by immunization of CCDC65 recombinant protein (Accession Number: NM_033124). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:

WB: 1:500-1:5000


IP: 1:500-1:5000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108984(CCDC65 Antibody) at dilution of

1:1000



IP Result of anti-CCDC65 (IP:Catalog No:108984, 4ug; Detection:Catalog No:108984 1:1000) with mouse testis tissue lysate 4000ug.