

# CCDC65 antibody



Catalog Number: 108984

## Product name

CCDC65 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CCDC65 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse testis tissue, rat testis tissue. Observed molecular weight by Western-blot: 57 kDa

## Preparation

This antibody was obtained by immunization of CCDC65 recombinant protein (Accession Number: NM\_033124). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:108984(CCDC65 Antibody) at dilution of



IP Result of anti-CCDC65 (IP:Catalog No:108984, 4ug; Detection:Catalog No:108984 1:1000) with mouse testis tissue lysate 4000ug.