

# CCDC67 antibody



Catalog Number: 108985

## Product name

CCDC67 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CCDC67 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 60-70 kDa

## Preparation

This antibody was obtained by immunization of CCDC67 recombinant protein (Accession Number: BC031247). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:108985(CCDC67 Antibody) at dilution of 1:600