

CCL25/TECK antibody



Catalog Number: 109002

Product name

CCL25/TECK antibody

Specificity

Human; other species not tested.

Antibody description

CCL25/TECK Rabbit Polyclonal antibody. Positive WB detected in Transfected HEK-293 cells.

Preparation

This antibody was obtained by immunization of CCL25/TECK recombinant protein (Accession Number: XM_017027121). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109002(CCL25 Antibody) at dilution of 1:700