

# CD23 antibody

Catalog Number: 109011

## Product name

CD23 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD23 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue.

## Preparation

This antibody was obtained by immunization of CD23 recombinant protein (Accession Number: BC064417). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:50-1:500

## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109011 (CD23, FCER2 Antibody) at dilution of 1:200 (under 10x lens). Heat-mediated antigen retrieval with Tris-EDTA buffer (pH 9).



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109011 (CD23, FCER2 Antibody) at dilution of 1:200 (under 40x lens). Heat-mediated antigen retrieval with Tris-EDTA buffer (pH 9).