

CD24 antibody



Catalog Number: 109013

Product name

CD24 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CD24 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HeLa cells, rat spleen tissue. Observed molecular weight by Western-blot: 30-70 kDa

Preparation

This antibody was obtained by immunization of CD24 recombinant protein (Accession Number: BC007674). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109013(CD24 antibody) at dilution of 1:1000