

CD244 antibody



Catalog Number: 109014

Product name

CD244 antibody

ELISA, IHC

Specificity

Human; other species not tested.

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Antibody description

CD244 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human small intestine tissue, human testis tissue.

Preparation

This antibody was obtained by immunization of CD244 recombinant protein (Accession Number: NM_016382). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality


Polyclonal

Ig Type


Rabbit IgG

Applications

Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109014(CD244 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109014(CD244 Antibody) at dilution of 1:50 (under 40x lens)