

# CD30 antibody



Catalog Number: 109022

## Product name

CD30 antibody

Recommended Dilution:

IHC: 1:50-1:500

## Specificity

Human; other species not tested.

## Antibody description

CD30 Rabbit Polyclonal antibody. Positive IHC detected in human lymphoma tissue, human tonsillitis tissue.

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_001243).  
Purification method: Antigen affinity purified.

## Formulation

0.1M NaHCO<sub>3</sub>, 0.1M glycine, 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


ELISA, IHC

## Dilutions

## Validations



Immunohistochemistry of paraffin-embedded human T-cell anaplastic large cell lymphoma (ALCL) tissue slide using Catalog No:109022(TNFRSF8 Antibody) at dilution of 1:200 (under 10x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).



Immunohistochemistry of paraffin-embedded human T-cell anaplastic large cell lymphoma (ALCL) tissue slide using Catalog No:109022(TNFRSF8 Antibody) at dilution of 1:200 (under 40x lens). Heat mediated antigen retrieved with Tris-EDTA buffer(pH9).