

# CD109 antibody

Catalog Number: 109040

## Product name

CD109 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD109 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells. Positive IP detected in PC-3 cells. Observed molecular weight by Western-blot: 180 kDa, 150 kDa, 120 kDa, 75 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: XM\_005248659). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109040(CD109 antibody) at dilution of 1:500



IP Result of anti-CD109 (IP:Catalog No:109040, 4ug; Detection:Catalog No:109040 1:600) with PC-3 cells lysate 5000ug.