

# CD11B/Integrin alpha M antibody



Catalog Number: 109041

## Product name

CD11B/Integrin alpha M antibody

## Specificity

Human; other species not tested.

## Antibody description

CD11B/Integrin alpha M Rabbit Polyclonal antibody. Positive FC detected in HL-60 cells. Positive IHC detected in human spleen tissue, human tonsillitis tissue. Positive WB detected in HL-60 cells, Jurkat cells. Observed molecular weight by Western-blot: 110-170 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_000632). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB, FC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109041(CD11B/Integrin alpha M Antibody) at dilution of 1:600



Immunohistochemical of paraffin-embedded human spleen using Catalog No:109041(CD11B antibody) at dilution of 1:50 (under 10x lens)




Immunohistochemical of paraffin-embedded human spleen using Catalog No:109041(CD11B antibody) at dilution of 1:50 (under 40x lens)

# CD11B/Integrin alpha M antibody



Catalog Number: 109041

---



1X10<sup>6</sup> HL-60 cells were stained with 0.2ug CD11B/Integrin alpha M antibody (Catalog No:109041, red) and control antibody (blue). Fixed with 4% PFA blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.