

CD11d/Integrin alpha D antibody



Catalog Number: 109044

Product name

CD11d/Integrin alpha D antibody

Specificity

Human; other species not tested.

Antibody description

CD11d/Integrin alpha D Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 126-150 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_005353). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109044(CD11d antibody) at dilution of 1:800