

CD133-1,2,3,5,7 antibody



Catalog Number: 109050

Product name

CD133-1,2,3,5,7 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD133-1,2,3,5,7 Rabbit Polyclonal antibody. Positive IP detected in Y79 cells. Positive WB detected in Y79 cells, mouse lung tissue. Observed molecular weight by Western-blot: 110 kDa

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_001145851). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP

Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



Y79 cells were subjected to SDS PAGE followed by western blot with Catalog No:109050(PROM1-1,2,3,5,7 antibody) at dilution



IP Result of anti-CD133-1,2,3,5,7 (IP:Catalog No:109050, 4ug; Detection:Catalog No:109050 1:300) with Y79 cells lysate 2000ug.