

# CD151 antibody

Catalog Number: 109056

## Product name

CD151 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CD151 Rabbit Polyclonal antibody. Positive IP detected in mouse spleen tissue. Positive WB detected in mouse spleen tissue. Observed molecular weight by Western-blot: 32 kDa

## Preparation

This antibody was obtained by immunization of CD151 recombinant protein (Accession Number: NM\_139029). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

## Validations



mouse spleen tissue were subjected to SDS PAGE followed by western blot with Catalog No:109056(CD151 antibody) at dilution of 1:300



IP Result of anti-CD151 (IP:Catalog No:109056, 3ug; Detection:Catalog No:109056 1:300) with mouse spleen tissue lysate 4000ug.