

# CCR3 antibody



Catalog Number: 109076

## Product name

CCR3 antibody

## Specificity

Human; other species not tested.

## Antibody description

CCR3 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue.

## Preparation

This antibody was obtained by immunization of CCR3 recombinant protein (Accession Number: NM\_001837). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, IHC

## Dilutions


Recommended Dilution:

IHC: 1:100-1:400

## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109076(CCR3 Antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109076(CCR3 Antibody) at dilution of 1:200 (under 40x lens)