

CCT6B antibody



Catalog Number: 109086

Product name

CCT6B antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CCT6B Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue. Observed molecular weight by Western-blot: 58kd

Preparation

This antibody was obtained by immunization of CCT6B recombinant protein (Accession Number: NM_006584). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:109086(CCT6B antibody) at dilution of 1:500