

# CCT6B antibody



Catalog Number: 109086

## Product name

Polyclonal

CCT6B antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

CCT6B Rabbit Polyclonal antibody. Positive WB detected in mouse ovary tissue. Observed molecular weight by Western-blot: 58kd

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:500-1:5000

This antibody was obtained by immunization of CCT6B recombinant protein (Accession Number: NM\_006584). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:109086(CCT6B antibody) at dilution of 1:500