

CDC14B antibody



Catalog Number: 109092

Product name

Polyclonal

CDC14B antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

CDC14B Rabbit Polyclonal antibody. Positive WB detected in SMMC-7721 cells. Observed molecular weight by Western-blot: 56 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:1000

This antibody was obtained by immunization of CDC14B recombinant protein (Accession Number: NM_033331). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



SMMC-7721 cells were subjected to SDS PAGE followed by western blot with Catalog No:109092(CDC14B Antibody) at dilution of 1:300