

# CDC25C antibody

Catalog Number: 109098

## Product name

CDC25C antibody

## Specificity

Human; other species not tested.

## Antibody description

CDC25C Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, HepG2 cells, Raji cells. Positive IP detected in K-562 cells. Observed molecular weight by Western-blot: 53 kDa

## Preparation

This antibody was obtained by immunization of CDC25C recombinant protein (Accession Number: XM\_047417952). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:500-1:5000

## Validations



K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:109098(CDC25C antibody) at dilution of 1:800



IP Result of anti-CDC25C (IP:Catalog No:109098, 4ug; Detection:Catalog No:109098 1:1000) with K-562 cells lysate 3600ug.