

# CDCA4 antibody



Catalog Number: 109115

---

## Product name

CDCA4 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CDCA4 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells. Observed molecular weight by Western-blot: 33kd

## Preparation

This antibody was obtained by immunization of CDCA4 recombinant protein (Accession Number: NM\_017955). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109115(CDCA4 antibody) at dilution of 1:500