

CD72 antibody



Catalog Number: 109137

Product name

CD72 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD72 Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue. Observed molecular weight by Western-blot: 40kd

Preparation

This antibody was obtained by immunization of CD72 recombinant protein (Accession Number: NM_001782). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109137(CD72 antibody) at dilution of 1:300