

CD80 antibody

Catalog Number: 109139

Product name

CD80 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CD80 Rabbit Polyclonal antibody. Positive FC detected in Raji cells. Positive WB detected in Human peripheral blood leukocyte cells, mouse skeletal muscle tissue, mouse thymus tissue, RAW 264.7 cells. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of CD80 recombinant protein (Accession Number: NM_005191). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



Human peripheral blood leukocyte cells were subjected to SDS PAGE followed by western blot with Catalog No:109139(B7-1 Antibody) at dilution of 1:300



1X10⁶ Raji cells were stained with 0.2ug B7-1 antibody (Catalog No:109139, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1000.