

CD81 antibody

Catalog Number: 109140

Product name

CD81 antibody

Specificity

Human; other species not tested.

Antibody description

CD81 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue.

Preparation

This antibody was obtained by immunization of CD81 recombinant protein (Accession Number: NM_004356). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:50-1:500

Validations



Immunohistochemical of paraffin-embedded human tonsillitis using Catalog No:109140(CD81 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human tonsillitis using Catalog No:109140(CD81 antibody) at dilution of 1:100 (under 40x lens)