

# CD84 antibody



Catalog Number: 109142

## Product name

CD84 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

CD84 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, mouse liver tissue, mouse spleen tissue, mouse thymus tissue, U-937 cells. Observed molecular weight by Western-blot: 45 kDa, 66-70 kDa

## Preparation

This antibody was obtained by immunization of CD84 recombinant protein (Accession Number: NM\_003874). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:109142(CD84 antibody) at dilution of 1:3000