

# CD89 antibody

Catalog Number: 109144

---

## Product name

CD89 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD89 Rabbit Polyclonal antibody. Positive WB detected in HL-60 cells. Observed molecular weight by Western-blot: 55-75 kDa

## Preparation

This antibody was obtained by immunization of CD89 recombinant protein (Accession Number: NM\_002000). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:1000

## Validations



HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:109144(FCAR antibody) at dilution of 1:300