

CD8b antibody



Catalog Number: 109146

Product name

Polyclonal

CD8b antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

CD8b Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HepG2 cells. Observed molecular weight by Western-blot: 37-43 kDa

Dilutions

Recommended Dilution:

Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of CD8b recombinant protein (Accession Number: NM_172213). Purification method: Antigen affinity purified.

Validations


Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109146(CD8B antibody) at dilution of 1:300