

# CD94 antibody



Catalog Number: 109149

## Product name

CD94 antibody

## Specificity

Human; other species not tested.

## Antibody description

CD94 Rabbit Polyclonal antibody. Positive IHC detected in human tonsillitis tissue. Positive WB detected in Jurkat cells, HL-60 cells. Observed molecular weight by Western-blot: 25-30 kDa

## Preparation

This antibody was obtained by immunization of CD94 recombinant protein (Accession Number: NM\_002262). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


## Validations



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109149(KLRD1 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:109149(KLRD1 Antibody) at dilution of 1:50 (under 40x lens)



(under 40x lens)

Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109149(KLRD1 antibody) at dilution of 1:500