

CDCP1 antibody



Catalog Number: 109151

Product name

CDCP1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CDCP1 Rabbit Polyclonal antibody. Positive WB detected in A549 cells. Observed molecular weight by Western-blot: 130 kDa

Preparation

This antibody was obtained by immunization of CDCP1 recombinant protein (Accession Number: NM_178181). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



A549 cells were subjected to SDS PAGE followed by western blot with Catalog No:109151(CDCP1 Antibody) at dilution of 1:600