

# CDH10-specific antibody



Catalog Number: 109152

## Product name

CDH10-specific antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CDH10-specific Rabbit Polyclonal antibody. Positive WB detected in rat brain tissue, mouse lung tissue, rat lung tissue. Observed molecular weight by Western-blot: 120 kDa

## Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM\_006727). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



rat brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:109152(CDH10-specific Antibody) at dilution of 1:1000