

CDX1 antibody



Catalog Number: 109176

Product name

CDX1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CDX1 Rabbit Polyclonal antibody. Positive WB detected in mouse colon tissue, mouse small intestine tissue, mouse stomach tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of CDX1 recombinant protein (Accession Number: BC096252). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse colon tissue were subjected to SDS PAGE followed by western blot with Catalog No:109176(CDX1 Antibody) at dilution of 1:600