

CDX4 antibody



Catalog Number: 109177

Product name

CDX4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CDX4 Rabbit Polyclonal antibody. Positive WB detected in SKOV-3 cells, mouse thymus tissue, Sp2/0 cells. Observed molecular weight by Western-blot: 37 kDa

Preparation

This antibody was obtained by immunization of CDX4 recombinant protein (Accession Number: NM_005193). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



SKOV-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:109177(CDX4 antibody) at dilution of 1:1000