

# CDYL antibody



Catalog Number: 109178

## Product name

CDYL antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CDYL Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, MCF7 cells, PC-3 cells. Observed molecular weight by Western-blot: 66 kDa, 60 kDa

## Preparation

This antibody was obtained by immunization of CDYL recombinant protein (Accession Number: NM\_004824). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

## Validations



mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109178(CDYL antibody) at dilution of 1:500