

CDYL antibody



Catalog Number: 109178

Product name

Polyclonal

CDYL antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, WB

CDYL Rabbit Polyclonal antibody. Positive WB detected in mouse testis tissue, MCF7 cells, PC-3 cells. Observed molecular weight by Western-blot: 66 kDa, 60 kDa

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000

Preparation

This antibody was obtained by immunization of CDYL recombinant protein (Accession Number: NM_004824). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

mouse testis tissue were subjected to SDS PAGE followed by western blot with Catalog No:109178(CDYL antibody) at dilution of 1:500