

CEPT1 antibody



Catalog Number: 109189

Product name

CEPT1 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CEPT1 Rabbit Polyclonal antibody. Positive IF detected in HepG2. Positive IHC detected in human pancreas.

Preparation

This antibody was obtained by immunization of CEPT1 recombinant protein (Accession Number: NM_001007794). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, IF


Dilutions

Recommended Dilution:


IHC: 1:20-1:200

IF: 1:10-1:100


Validations



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:109189(CEPT1 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas using Catalog No:109189(CEPT1 antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HepG2 cells, using CEPT1 antibody Catalog No:109189 at 1:25 dilution and FITC- labeled goat anti-Rabbit IgG (green).