

# CERCAM antibody



Catalog Number: 109190

## Product name

Polyclonal

CERCAM antibody

## Ig Type

## Specificity

Rabbit IgG

human, mouse, rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

CERCAM Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 68kd,48kd

## Dilutions


Recommended Dilution:

WB: 1:500-1:5000

## Preparation

This antibody was obtained by immunization of CERCAM recombinant protein (Accession Number: BC011811). Purification method: Antigen affinity purified.

## Validations



## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109190(CERCAM antibody) at dilution of 1:1500

## Clonality