

CES1 antibody



Catalog Number: 109192

Product name

CES1 antibody

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

CES1 Rabbit Polyclonal antibody. Positive IHC detected in human liver tissue, human kidney tissue, human ovary tissue. Positive WB detected in HEK-293 cells, Jurkat cells. Observed molecular weight by Western-blot: 62kd




Preparation

This antibody was obtained by immunization of CES1 recombinant protein (Accession Number: NM_001266). Purification method: Antigen affinity purified.

HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109192(CES1 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT


Clonality

Polyclonal

Immunohistochemical of paraffin-embedded human liver using Catalog No:109192(CES1 antibody) at dilution of 1:200 (under 10x lens)

Ig Type

Rabbit IgG



Applications

ELISA, WB, IHC

Immunohistochemical of paraffin-embedded human liver using Catalog No:109192(CES1 antibody) at dilution of 1:200 (under 40x lens)

Dilutions

Recommended Dilution:

WB: 1:200-1:2000