

CES2 antibody



Catalog Number: 109193

Product name

CES2 antibody

Specificity

Human; other species not tested.

Antibody description

CES2 Rabbit Polyclonal antibody. Positive WB detected in human stomach tissue, human liver tissue. Positive IHC detected in human colon cancer tissue. Observed molecular weight by Western-blot: 60kd

Preparation

This antibody was obtained by immunization of CES2 recombinant protein (Accession Number: BC032095). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



human stomach tissue were subjected to SDS PAGE followed by western blot with Catalog No:109193(CES2 antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109193(CES2 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109193(CES2 antibody) at dilution of 1:100 (under 40x lens)