

CES3 antibody

Catalog Number: 109194

Product name

CES3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

CES3 Rabbit Polyclonal antibody. Positive IHC detected in human colon cancer tissue. Positive WB detected in HeLa cells, human liver tissue. Observed molecular weight by Western-blot: 60-62kd

Preparation

This antibody was obtained by immunization of CES3 recombinant protein (Accession Number: NM_001185177). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109194(CES3 antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human colon cancer using Catalog No:109194(CES3 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human colon cancer using Catalog

CES3 antibody



Catalog Number: 109194

No:109194(CES3 antibody) at dilution of 1:100 (under 40x lens)