

# CFP antibody

Catalog Number: 109200

## Product name

CFP antibody

## Specificity

Human; other species not tested.

## Antibody description

CFP Rabbit Polyclonal antibody. Positive IHC detected in human spleen tissue. Positive WB detected in Jurkat cells. Observed molecular weight by Western-blot: 50 kDa

## Preparation

This antibody was obtained by immunization of CFP recombinant protein (Accession Number: NM\_002621). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

## Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109200(CFP Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human spleen tissue slide using Catalog No:109200(CFP Antibody) at dilution of 1:50 (under 40x lens)