

CGRRF1 antibody



Catalog Number: 109204

Product name

CGRRF1 antibody

Specificity

Human, Rat; other species not tested.

Antibody description

CGRRF1 Rabbit Polyclonal antibody. Positive WB detected in rat kidney tissue, L02 cells. Observed molecular weight by Western-blot: 42-48 kDa

Preparation

This antibody was obtained by immunization of CGRRF1 recombinant protein (Accession Number: NM_006568). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



rat kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109204(CGRRF1 Antibody) at dilution of 1:600