

# CGRRF1 antibody



Catalog Number: 109204

## Product name

Polyclonal

CGRRF1 antibody

## Ig Type

## Specificity

Rabbit IgG

Human, Rat; other species not tested.

## Applications

## Antibody description

ELISA, WB

CGRRF1 Rabbit Polyclonal antibody. Positive WB detected in rat kidney tissue, L02 cells. Observed molecular weight by Western-blot: 42-48 kDa

## Dilutions

Recommended Dilution:

## Preparation

WB: 1:200-1:2000

This antibody was obtained by immunization of CGRRF1 recombinant protein (Accession Number: NM\_006568). Purification method: Antigen affinity purified.

## Validations


## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality



rat kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:109204(CGRRF1 Antibody) at dilution of 1:600