

CHAC2 antibody



Catalog Number: 109205

Product name

CHAC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHAC2 Rabbit Polyclonal antibody. Positive IF detected in MCF-7 cells. Positive WB detected in Jurkat cells, K-562 cells, mouse brain tissue, Raji cells. Observed molecular weight by Western-blot: 21kd

Preparation

This antibody was obtained by immunization of CHAC2 recombinant protein (Accession Number: NM_001008708). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109205(CHAC2 antibody) at dilution of 1:500



Immunofluorescent analysis of MCF-7 cells, using CHAC2 antibody Catalog No:109205 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).