

CHAD antibody



Catalog Number: 109206

Product name

CHAD antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHAD Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 41 kDa

Preparation

This antibody was obtained by immunization of CHAD recombinant protein (Accession Number: NM_001267). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:109206(CHAD antibody) at dilution of 1:400