

CHAF1A antibody



Catalog Number: 109207

Product name

CHAF1A antibody

IP: 1:500-1:5000

IHC: 1:20-1:200

Specificity

Human, Mouse; other species not tested.

Antibody description

CHAF1A Rabbit Polyclonal antibody. Positive IHC detected in human cervix tissue. Positive WB detected in Jurkat cells, HeLa cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 150 kDa

Preparation

This antibody was obtained by immunization of CHAF1A recombinant protein (Accession Number: NM_005483). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:1000-1:10000


Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:109207(CHAF1A antibody) at dilution of 1:2000



IP Result of anti-CHAF1A (IP:Catalog No:109207, 4ug; Detection:Catalog No:109207 1:1000) with Jurkat cells lysate 2400ug.



Immunohistochemistry of paraffin-embedded human cervix tissue slide using Catalog

CHAF1A antibody



Catalog Number: 109207

No:109207(CHAF1A Antibody) at dilution of 1:50 (under 40x lens)