

CHCHD1 antibody



Catalog Number: 109210

Product name

CHCHD1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHCHD1 Rabbit Polyclonal antibody. Positive IF detected in A431 cells. Positive WB detected in A431 cells. Observed molecular weight by Western-blot: 13 kDa

Preparation

This antibody was obtained by immunization of CHCHD1 recombinant protein (Accession Number: NM_203298). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IF: 1:50-1:500

Validations



A431 cells were subjected to SDS PAGE followed by western blot with Catalog No:109210(CHCHD1 antibody) at dilution of 1:300



Immunofluorescent analysis of (10% Formaldehyde) fixed A431 cells using Catalog No:109210(CHCHD1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)