

CHCHD6 antibody



Catalog Number: 109215

Product name

CHCHD6 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200


Specificity

Human; other species not tested.

Validations

Antibody description

CHCHD6 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Positive IHC detected in human stomach tissue. Observed molecular weight by Western-blot: 26-29 kDa




Preparation

This antibody was obtained by immunization of CHCHD6 recombinant protein (Accession Number: NM_032343). Purification method: Antigen affinity purified.

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:109215(CHCHD6 antibody) at dilution of 1:500

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.



Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Immunohistochemical of paraffin-embedded human stomach using Catalog No:109215(CHCHD6 antibody) at dilution of 1:100 (under 40x lens)

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, IP



Dilutions

Recommended Dilution:

WB: 1:1000-1:10000

Immunohistochemical of paraffin-embedded

CHCHD6 antibody

Catalog Number: 109215

human stomach using Catalog
No:109215(CHCHD6 antibody) at dilution of 1:100

IP Result of anti-CHCHD6 (IP:Catalog No:109215,
4ug; Detection:Catalog No:109215 1:500) with
HeLa cells lysate 1200ug.

