

CHDH antibody



Catalog Number: 109222

Product name

CHDH antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHDH Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, HepG2 cells, mouse kidney tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 60 kDa

Preparation

This antibody was obtained by immunization of CHDH recombinant protein (Accession Number: NM_018397). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IHC: 1:20-1:200


Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109222(CHDH antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human liver using Catalog No:109222(CHDH antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:109222(CHDH antibody) at dilution of 1:100 (under 40x lens)