

CHI3L1 antibody



Catalog Number: 109225

Product name

Polyclonal

CHI3L1 antibody

Ig Type

Specificity

Rabbit IgG

Human, Mouse, Rat; other species not tested.

Applications

Antibody description

ELISA, IHC

CHI3L1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue, human liver tissue.

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Preparation

This antibody was obtained by immunization of CHI3L1 recombinant protein (Accession Number: NM_001276). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:109225(CHI3L1 antibody) at dilution of 1:50 (under 40x lens)