

CHI3L1 antibody

Catalog Number: 109225

Product name

CHI3L1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

CHI3L1 Rabbit Polyclonal antibody. Positive IHC detected in human ovary tumor tissue, human liver tissue.

Preparation

This antibody was obtained by immunization of CHI3L1 recombinant protein (Accession Number: NM_001276). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human ovary tumor using Catalog No:109225(CHI3L1 antibody) at dilution of 1:50 (under 40x lens)