

# CHMP7 antibody



Catalog Number: 109237

## Product name

CHMP7 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

CHMP7 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue. Positive WB detected in HEK-293 cells. Observed molecular weight by Western-blot: 51-55kd

## Preparation

This antibody was obtained by immunization of CHMP7 recombinant protein (Accession Number: NM\_152272). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, IHC, WB


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109237(CHMP7 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:109237(CHMP7 antibody) at dilution of 1:100 (under 40x lens)



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:109237(CHMP7 Antibody) at dilution of 1:1000