

CHN2 antibody



Catalog Number: 109239

Product name

CHN2 antibody

Specificity

human, mouse, rat; other species not tested.

Antibody description

CHN2 Rabbit Polyclonal antibody. Positive WB detected in mouse cerebellum tissue. Observed molecular weight by Western-blot: 46-50kd

Preparation

This antibody was obtained by immunization of CHN2 recombinant protein (Accession Number: NM_004067). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:109239(CHN2 antibody) at dilution of 1:600